



Skedco HMH SKED® Strap Kit

Patient Care | Flexible Litters and Accessories



The HMH Sked® is designed for use in mass casualty incidents or Haz-Mat incidents. Patients can be rolled or dragged onto the HMH Sked® rather than lifted. This creates less rescuer fatigue and requires fewer rescuers to be exposed in the Hot Zone at a Haz-Mat incident. The HMH Sked® can be easily dragged by personnel wearing level A protection or even pulled by a rope connected outside the Hot Zone.

Key Features: No grommets to trap contaminants • Stores flat in stacks of up to 200 • Ideal for use in mass casualty incidents • Polypropylene straps are chemical resistant • 2" Side release buckles • Not for lifting patients

Additional Information

Shipping weight | 2 lbs

Product Code

SK252

