



Included are 4 HMH Skeds rolled and stored in rapid deployment bags. Also included are 8 Shuttle Sked rope kits in a single rapid deployment bag to keep everything organized. It would be wise to have up to 4 more HMD Skeds to move patients from the contaminated Shuttle Skeds to decon. HMH Shuttle Sked ropes are constructed of 1/2 inch solid braid polypropylene. They are 200ft. long with a figure 8 knot saftied with a barrel knot and taped to prevent unraveling while in use. A Large aluminum D carabiner is attached. The rope kit includes a Haz-Mat rope bag for quick tangle free deployment. These ropes are attached to both ends of a HMD Sked for dragging into and out of the Hot-zone.

Additional Information

Shipping weight | 132 lbs Brand | Skedco

Product Code

SK2510

