



HMH SKED® Rescue System w/strap kit Patient Care | Flexible Litters and Accessories



The HMH Sked® is designed for use in mass casualty incidents or Haz-Mat incidents. Patients can be rolled or dragged onto the HMH Sked® rather than lifted. This creates less rescuer fatigue and requires fewer rescuers to be exposed in the Hot Zone at a Haz-Mat incident. The HMH Sked® can be easily dragged by personnel wearing level A protection or even pulled by a rope connected outside the Hot Zone.

Key Features

- No grommets to trap contaminants
- Stores flat in stacks of up to 200
- Ideal for use in mass casualty incidents
- Polypropylene straps are chemical resistant
- 2" Side release buckles
- Not for lifting patients

Additional Information

Shipping weight | 13.3 lbs

Weight | 3.00 lbs (5.900 kg)

Compliant | CE

Tags

[Made in USA](#)

Product Code

SK250

