## HALF SKED® PATIENT DRAG



Made from the same durable Sked® material, but only as half as long as the original Sked®. It allows patients to bend at the hips, providing greater flexibility in the most demanding and tight spaces. Not to be used as a hoisting or immobilization device.

## **Key Features:**

Perfect for tight spaces • Allows greater flexibility • 2" wide chest and leg straps with easy release buckles

Weight: 8.00 lbs (3.700 kg)

Orange **SK220** 

CERTIFIED CE



Pigeon Mountain Industries, Inc. | PO Box 803, LaFayette, GA 30278 | Fax: 706-764-1531 | Email: custserv@pmirope.com



## **NOTES**



Pigeon Mountain Industries, Inc. | PO Box 803, LaFayette, GA 30278 | Fax: 706-764-1531 | Email: custserv@pmirope.com



